



MULTINEB SPECIFICATIONS

Sample Connection: Direct

Technology: Flow-blurring / Neutral Suction

Self Aspirating: No (No back-pressure)

Gas Connector: Quick release

Material: PFA, ETFE, PEEK and LCP

UNSPSC Code: 41113000

Chemical Compatibility: Aqueous, HF acid and Organic solvents

Notes: Ideal for routine analysis of samples up to 20 % TDS with good precision and better sensitivity than a SeaSpray; inert construction makes it suitable for most solution types.

TDS Handling Capability: 20 %TDS

Flow rate: 0.25-1ml/min (for each capillary)

Pressure range: 180-200 Kpa (0.7-1 lpm Ar)

Other Features: Not suitable for use with SVS/AVS