



# ONENEB 2 SPECIFICATIONS

---

---

**Sample Connection:** Direct (1/16" PFA tubing)

**Technology:** Flow-blurring

**Self Aspirating:** No

**Gas Connector:** Quick release (4mm)

**Material:** PFA, ETFE, PEEK and LCP

**UNSPSC Code:** 41113000

**Chemical Compatibility:** Aqueous, HF acid and Organic solvents

**Notes:** Ideal for routine analysis of samples up to 25 % TDS with good precision and better sensitivity than a SeaSpray; inert construction makes it suitable for virtually any solution. Delivers outstanding precision: typically < 1% RSD.

**TDS Handling Capability:** 25 %TDS

**Flow rate:** 0.5-2ml/min (1ml/min recommended)

**Pressure range:** 150 Kpa (0.7 lpm Ar)

**Other Features:** Not suitable for use with SVS/AVS