



Oneneb 2 specifications

Sample Connection: Direct (1/16" PFA tubing)

Technology: Flow-blurring

Self Aspirating: No

Gas Connector: Quick release (4mm)

Material: PFA, ETFE, PEEK and LCP

UNSPSC Code: 41113000

Chemical Compatibility: Aqueous, HF acid and Organic solvents

Notes: Ideal for routine analysis of samples up to 25 % TDS with good precision and better sensitivity than a SeaSpray; inert construction makes it suitable for virtually any solution. Delivers outstanding precision: typically < 1% RSD.

TDS Handling Capability: 25 %TDS

Flow rate: 0.5-2ml/min (1ml/min recommended)

Pressure range: 150 Kpa (0.7 lpm Ar)

Other Features: Not suitable for use with SVS/AVS